

Instruction manual Blend Tag® and Blend Tag® -Plus- 0810

F0938K

Blend Taq[®] / Blend Taq[®] -Plus-

<Blend Taq[®]> <Blend Taq[®] -Plus->

BTQ-101 250 U 200 reactions BTQ-201 250 U 200 reactions

Store at -20°C

Contents

- [1] Introduction
- [2] Components
- [3] Quality testing
- [4] Primer design
- [5] Cloning of PCR products
- Protocol [6]
 - 1. Standard reaction setup
 - 2. Cycling conditions
- Examples [7]
- [8] Troubleshooting
- [9] Related products
- [10] References

CAUSION

All reagents in this kit are intended for research purposes. Do not use for diagnosis or clinical purposes. Please observe general laboratory precaution and utilize safety while using this kit.

- Blend Taq[®] is a registered trademark of Toyobo Co., Ltd. in Japan.

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp

CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



[1] Introduction

Description

Blend Taq[®] and Blend Taq[®] -Plus- are highly efficient Taq-based DNA polymerases developed based on the Barns' method.¹⁾ This method uses a DNA polymerase v lacked $3' \rightarrow 5'$ exonuclease (proofreading) activity (e.g., Taq DNA polymerase) a small amount of an archaeal DNA polymerase with proofreading activity. Because une proofreading activity repairs misincorporated nucleotide bases causing the termination of the polymerase reaction, PCR with a 'mixed' enzyme solution enables highly efficient amplification.

The enzyme solution of Blend Taq[®] -Plus- contains anti-Taq DNA polymerase antibodies that inhibit polymerase activity, allowing for Hot Start PCR.

Blend Taq[®] and Blend Taq[®] -Plus- generate dA overhang-ended PCR products. Therefore, the PCR products can be cloned using a standard TA-cloning method.

Features

-This enzyme is effective for the amplification of various targets from small template amounts. The elongation ability of this enzyme is much greater than that of the normal Taq DNA polymerase.

-Hot Start technology using anti-Taq DNA polymerase antibodies results in highly efficient amplification. <Blend Taq[®] (Code No. BTQ-101) does not use hot start>

-The PCR error ratio of this enzyme is approximately 3-4 times less than that of Taq DNA polymerase.




Fig. 1. The principles of the Barns' technology.

[2]	Components	The provided reagents include the following components for 200 reactions:		
	<blend taq<sup="">®></blend>	Blend Taq [®] (2.5 U/μl) 10x PCR Buffer for Blend Taq [®] 2 mM dNTPs	100 μl 1.0 ml 1.0 ml	
	<blend taq<sup="">® -Plus-></blend>	Blend Taq [®] -Plus- (2.5 U/µl)* 10x PCR Buffer for Blend Taq [®] 2 mM dNTPs *This enzyme solution contains anti-Tac	100 μl 1.0 ml 1.0 ml I DNA polymerase antibodies.	

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



[3]	Quality Testing	Quality check was performed by amplifying the human β -globin gene (17.5 kb).
[4]	Primer Design	PCR primers should be designed according to general guidelines. For the amplification of a long target (≥ 6 kb), the melting temperature (Tm) of the primers should be set over 70 °C.
[5]	Cloning of PCR products	The PCR products of Blend $Taq^{\ensuremath{\mathbb{R}}}$ and Blend $Taq^{\ensuremath{\mathbb{R}}}$ -Plus- can be cloned according to a standard TA-cloning method.

[6] **Protocol** 1. Standard reaction setup

The following procedures are designed to be used with the components provided in this kit. Before preparing the reaction mixture, all components should be completely thawed, except for the enzyme solution.

Component	Volume	Final Concentration
10x Buffer	5 µl	1 x
2 mM dNTPs*	5 µl	0.2 mM each
10 pmol/µl Primer #1	1 µl	0.2 μΜ
10 pmol/µl Primer #2	1 µl	0.2 μΜ
Template DNA	Xμl {	Genomic DNA 10-1000 ng/50 µl Plasmid DNA 1-50 ng/50 µl cDNA ≤200 ng (RNA equiv.)/50 µl - <i>E. coli</i> cells (small amount)
PCR grade water	Υµl	
Blend Taq [®] (2.5 U/ μ l) or	0.5 µl	1.25 U / 50 μl
Blend Taq [®] -Plus (2.5 U/µl)		
Total reaction volume	50 µl	

* Do not use dNTPs from other kits or companies.

Notes:

- -Thin-wall tubes are recommended for PCR use. A total reaction volume of 50 μl is recommended.
- -For amplification of a long target (\geq 10 kb), the final concentration of dNTPs should be 0.3-0.4 mM.

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp

CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



2. Cycling conditions

The following cycling steps are recommended.

(1) Cycling conditions for < 6kb targets.

	_	
94 °C, 2 min	_	
94 °C, 30 sec.	-	
Tm-[5-10] °C*, 30 sec.		25-35 cycles
72 °C, 1 min/kb		
	94 °C, 30 sec. Tm-[5-10] °C*, 30 sec.	94 °C, 30 sec. Tm-[5-10] °C*, 30 sec.

*Tm value of the primer minus 5°C-10°C

(2) Cycling conditions for \geq 6kb products.

< 2-step cycle >			
Pre-denaturation	94 °C, 2 min		
Denaturation	94 °C, 30 sec.	•	25.25 1
Extension	68 °C, 1 min/kb		25-35 cycles


Note: If the Tm value of the primer is under 73 °C, the 3-step cycle is recommended.

Notes:

-Extension time should be set at 1min per 1 kb of target length. -For colony-direct PCR, the pre-denaturation step should be set for 4 min.

[7] Examples

Example 1. Comparison of the sensitivity of PCR with a Taq-based PCR enzyme.



Template: Human Genomic DNA Target: β-globin 3.6 kb


Lane 1: 0 ng human genomic DNA Lane 2: 5 ng human genomic DNA Lane 3: 10 ng human genomic DNA Lane 4: 20 ng human genomic DNA Lane 5: 40 ng human genomic DNA

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp

CHINA

TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140





Example 2. Colony-direct PCR using E. coli cells.

Lane 1: Insert (+) Lane 2: Insert (+) Lane 3: Insert (-) Lane 4: Insert (+)

F primer: TCGAGGTCGACGGTATCGAT (20mer GC=55%) R primer: CGCTCTAGAACTAGTGGATC (20mer GC=50%)

*The primers are designed on the vector.

[8] Troubleshooting

Symptoms	Cause	Solution
	Cycling conditions are not optimal	Lower annealing temperature increments to a maximum of Tm-5 °C.
		Increase the number of cycles by 2-5 cycles.
No PCR product / low yield	Primers are not good	Check the design and/or quality of primers.
	Too much <i>E. coli</i> cells (Colony direct PCR)	Decrease the amount of <i>E. coli</i> cells from colonies.
Smearing / extra band(s)	Cycling conditions are not optimal	Decrease the number of cycles by 2-5 cycles.

[9] Related products

Product name	Package	Code No.
Highly efficient cDNA synthesis kit	100 rxns	FSK-101
ReverTra Ace [®] -α-		
Highly efficient reverse transcriptase	10,000 U	TRT-101
ReverTra Ace [®]		
RNase inhibitor (Recombinant type)	2,500 U	SIN-201
High Speed PCR enzyme	200 U x1	LDP-101
KOD Dash		
High fidelity PCR enzyme	200 U x1	KOD-201
KOD -Plus-		
Highly reliable PCR enzyme	200 U x1	KFX-101
KOD FX		

[10] References

1) Barns WM, Proc. Natl. Acad. Sci. USA, 91: 2216-2220 (1994)

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp

CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



NOTICE TO PURCHASER: LIMITED LICENSE

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as the patented 5' Nuclease Process claims in US Patents Nos. 5,210,015 and 5,487,972), no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp **CHINA** TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140