



T4 DNA Ligase

LGA-111 1 x 400 units

Content : T4 DNA ligase (1-5 U/ul)
10 x buffer*

* rATP is not supplied. Must be added separately to ligation reactions

- Application

DNA Fragment and Nick Joining

T4 DNA Ligase can be used in the insertion of target genes into vectors or in reactions such as the addition of linker to DNA.

- Description

T4 DNA Ligase catalyzes the formation of phosphodiester bonds between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA and requires Mg^{2+} and rATP as catalysts.

It works on blunt ends, protruding ends and nicked sites of stranded DNA, but does little joining of RNA-RNA or RNA-DNA. Also, the ligation reaction rate varies depending on terminal structure and DNA sequence.

- Principle

Ligation activity

- One-Point Advice

- Blunt-End Ligation

Blunt-end ligation is slower than protruding end ligation (K_m value is approximately 100-fold higher in the former). If ligating blunt ends, DNA concentration should be increased, and 2-5 times the amount of enzyme for protruding end ligation used.

- Reaction Temperature

This enzyme, with an optimum temperature of $37^{\circ}C$, has poor thermostability. Long reactions should take place at $16^{\circ}C$. 1-2 hour reactions are possible at room temperature.

- EDTA

Since this enzyme requires Mg^{2+} as a cofactor, EDTA may inhibit the reaction as it forms cholate with Mg^{2+} ions.

When using DNA dissolved in buffer with a high EDTA concentration as the sample for ligation, we recommend substitution with sterile water or TE Buffer.



■ Amount of DNA

The molar ratio of plasmid DNA to insert DNA should be 1:3. For ligation to cosmid or phage, we recommend adjustment of the vector DNA-insert DNA molar ratio to 1:1, along with an increase in DNA concentration (more than 0.05-0.1 µg/µl)

● Basic Conditions for Reaction:

■ Ligation into Vector

Vector DNA	50-500ng
Insert DNA	3X Vector DNA
10X Buffer	2µl
10mM rATP	2µl
T4 DNA Ligase	
Protruding End	0.5-1U
Blunt End	1-5U

Total Volume 20µl

16°C, 4-16hr.

■ Linker Ligation

DNA	0.1-1µg
Linker	10X DNA quantity
10X Buffer	2µl
10mM rATP	2µl
T4 DNA Ligase	1-5U

Total Volume 20µl

16°C, 8-16hr.

◇ The reaction varies depending on the concentration and structure, base sequences, and DNA length.